


ZSI 100 FTM

MALLEABLE PENILE IMPLANT



A solution to to provide erection
to phalloplasty

- Designed for FTM morphology.
- For more information, visit www.zsimplants.ch

ZSI 100 FTM MALLEABLE PENILE IMPLANT FOR PHALLOPLASTY


DISTAL PART
Adjustable by cutting 5 mm by 5 mm.

CENTRAL PART
Cylinder of 22 mm in width. Designed for comfort and pleasant appearance.

PROXIMAL PART
30 mm by 25 mm. Is fixed securely onto the pubic bone.


STANDARD VERSION

- 1 component Penile Implant.
- Sturdy fixation on pubic bone to provide excellent stability and cosmetic results.
- Proximal part with 4 attachments.
- 1 cylinder adjustable from 13 cm to 20 cm.




NARROW VERSION

- 1 component Penile Implant.
- sturdy fixation on pubic bone to provide excellent stability and cosmetic results.
- Proximal part with 4 attachments.
- 1 cylinder adjustable from 13 cm to 22 cm.
- Diameter: 13mm.





SCALE 1:1



Proximal part made of stainless steel and silicone fixed securely onto the pubic bone with 4 attachments

It fits securely to the pubic bone



Optimal comfort while wearing trousers

