


ZSI 475 FTM

INFLATABLE PENILE IMPLANT



A solution to bring erection
to Phalloplasty

- Designed for FTM morphology.
- For more information, visit www.zsimplants.ch

ZSI 475 FTM INFLATABLE PENILE IMPLANT

DESIGNED TO PROVIDE EXCELLENT STABILITY AND COSMETIC RESULTS

ONE CYLINDER, 4 SIZES, 120MM, 150MM, 180MM AND 210MM IN LENGTH
EASY TO IMPLANT AND STURDY FIXATION

PROXIMAL PART
30mm by 25mm

CENTRAL PART
main body diameter \varnothing 21-22 mm

DISTAL PART
glans diameter 25 mm




Deflated cylinder



Realistically shaped glans



Reservoir implanted in the pelvic cavity




Flexible proximal part made of stainless steel and silicone fixed securely onto the pubic bone




Inflatable testis prosthesis integrated in the pump unit


TO FUNCTION




1. Press the pump unit several times to inflate the prosthesis.



2. Penis erected.



3. Keep pressed the deflating valves until cylinder is empty.



4. Penis flaccid.