






For more information, visit
www.zsimplants.ch




The best solution
for sex reassignment
surgeries




ZSI 475 FTM INFLATABLE PENILE IMPLANT FOR PHALLOPLASTY

- 3 components Penile Implant.
- Sturdy fixation on pubic bone.
- 1 cylinder, 4 sizes, 120/150/180/210 mm.
- 1 cylinder designed to provide excellent stability and cosmetic result.
- 25 mm wide realistically shaped glands stopper. Urethra passage is preserved.

| | | |
|-----------------------|------------------------------|----------------------|
| DISTAL PART | CENTRAL PART | PROXIMAL PART |
| Glans diameter 25 mm. | Main body diameter 21-22 mm. | 30 mm by 25 mm. |




- Pump designed as testis prosthesis for better cosmetic results.
- Easy access Deflation Button. *




Proximal part made of stainless steel and silicone fixed securely onto the pubic bone with 8 attachments

It fits securely on the pubic bone



Best confort while wearing trousers




ZSI 100 FTM MALLEABLE PENILE IMPLANT FOR PHALLOPLASTY

STANDARD VERSION


- 1 component Penile Implant.
- Sturdy fixation on pubic bone to provide excellent stability and cosmetic results.
- Proximal part with 4 attachments.
- 1 cylinder adjustable from 13 cm to 20 cm.

| | | |
|-------------------------------------|---|--|
| DISTAL PART | CENTRAL PART | PROXIMAL PART |
| Adjustable by cutting 5 mm by 5 mm. | Cylinder of 22 mm in width. Designed for comfort and pleasant appearance. | 30 mm by 25 mm. Is fixed securely onto the pubic bone. |




NARROW VERSION

- 1 component Penile Implant.
- sturdy fixation on pubic bone to provide excellent stability and cosmetic results.
- Proximal part with 4 attachments.
- 1 cylinder adjustable from 13 cm to 22 cm.
- Diameter: 13mm.




SCALE 1:1




Proximal part made of stainless steel and silicone fixed securely onto the pubic bone with 4 attachments

It fits securely to the pubic bone




Optimal confort while wearing trousers



ZSI 100 D4 MALLEABLE PENILE IMPLANT FOR METOIDIOPLASTY

- A pair of rods made of silicone with silver cable.
- Length 12 cm that can be reduced to 5 cm.
- Diameter 4 mm.



ZSI 200 EXPANDER FOR VAGINOPLASTY

- Inflatable expander.
- Made of biocompatible silicone.
- Length 90mm and 120mm.
- Diameter inflated: 40mm.

