


**RETRACTOR + SET OF MALE STAYS
CORPORA DILATORS
FURLOW MEASURER
CORPORA T HANDLE MEASURER
DOUBLE BLADE CARVERNOTOMES**



RETRACTOR + SET OF MALE STAYS

Qty 10 sets per box, disposable

2 sizes of stays:

- Blunt stays of 5mm: Réf. JUNE3000B5
- blunt stays of 12mm: Réf. JUNE 3000B12




8 blunt stays
(=8 Units in 1 Set)

+



1 Retractor Ring




CORPORA DILATOR

Qty 1 set: Ø 8 mm, 9 mm, 10 mm, 11 mm, 12 mm, 13 mm, 14 mm
Reusable, reference G-04



FURLOW MEASURER HEELAY

One per box
Reusable, reference H11640R



CORPORA T HANDLE MEASURER

Qty 1- Ø 8 mm - Blunt
Reusable, reference G-03



DOUBLE BLADE CAVERNOTOMES

Qty 1 set: Ø 6 mm, 7 mm, 8 mm, 9 mm, 10 mm, 11 mm, 12 mm, 14 mm per box
Reusable, reference G-01

